

Column: **MEGA-PAH FAST - 0.10mm, 0.10µm, 8m**
Catalog Code: F-PAH-010-010-8

Conditions:

Injection: Split 250°C, 400mL/min split flow,
0.50µL injection volume,
0.25µL injection volume,
0.10µL injection volume.
Detector: FID 340°C.
Oven Program: 80°C, 15°C/min, 330°C.
Carrier Gas: Hydrogen, 180kPa.

Sample:

8270 Calibration Mix #5, Revised 2000µg/mL in
Methylene Chloride.


- | | |
|------------------------|----------------------------|
| 1. Naphthalene | 10. Pyrene |
| 2. 1-methylnaphthalene | 11. Benzo[a]anthracene |
| 3. 2-methylnaphthalene | 12. Chrysene |
| 4. Acenaphthylene | 13. Benzo[b]fluoranthene |
| 5. Acenaphthene | 14. Benzo[k]fluoranthene |
| 6. Fluorene | 15. Benzo[a]pyrene |
| 7. Phenanthrene | 16. Indeno[1,2,3-cd]pyrene |
| 8. Anthracene | 17. Dibenzo[a,h]anthracene |
| 9. Fluoranthene | 18. Benzo[g,h,i]perylene |

PAHs - improve sensitivity with FAST-GC

FAST-GC
solutions



Column: **MEGA-PAH FAST - 0.10mm, 0.10µm, 8m**
Catalog Code: F-PAH-010-010-8



Sample:
8270 Calibration Mix #5, Revised 2000µg/mL in
Methylene Chloride.

0.82ng for each compound
into the column
0.33ng for each compound
into the column



**FAST-GC narrower
peaks allow you to
increase the sensitivity
level of your GC system**

MEGA®

www.mega.mi.it



since
1980



contact us: info@mega.mi.it

 follow us [@MEGAColumns](https://twitter.com/MEGAColumns)

www.mega.mi.it