

Column: **MEGA-FFAP EXT FAST** - 0.15mm, 0.25 μ m, 10m

Catalog Code: F-FFAPEXT-015-025-10

New Column

Conditions:

Injection: Split 250°C, 25mL/min split flow,
1.5 μ L injection volume.


Detector: FID 250°C.

Oven Program: 35°C (0.5min), 20°C/min.

Carrier Gas: Helium, 100kPa.

I. Ethanol

• = Ethanol impurities and other residual solvents



Column: **MEGA-FFAP EXT FAST** - 0.15mm, 0.25 μ m, 10m

Catalog Code: F-FFAPEXT-015-025-10

New Column

Conditions:

Injection: Split 250°C, 25mL/min split flow,
1.5 μ L injection volume.

Detector: FID 250°C.

Oven Program: 35°C (0.5min), 20°C/min.

Carrier Gas: Helium, 100kPa.

I. Ethyl Acetate

• = Ethyl Acetate impurities and other residual solvents

100ppm of Acetic Acid
in Ethyl Acetate



since
1980



contact us: info@mega.mi.it

 follow us [@MEGAColumns](https://twitter.com/MEGAColumns)

www.mega.mi.it