



2-methyl Butyric Acid

Chiral Separation on MEGA-DEX Columns

Column: **MEGA-DEX DMT Beta - 0.25mm, 0.25 μ m, 25m**
Catalog Code: DEX-DMT-B-025-025-25

Conditions:

Injection: Split 250°C, 0.5 μ L, 350 mL/min Split Flow.

Sample dilution: 10 μ L/1mL each in Methanol.

Detector: FID 250°C.

Oven Temp.: 50°C, 2°C/min, 150°C.

Carrier Gas: Hydrogen, 80kPa.

Contact us for additional information.

