



Alpha-Pinene, Linalool, Menthol (+/-) Chiral Separation on MEGA-DEX Columns

Column: **MEGA-DEX DET Beta - 0.25mm, 0.25 μ m, 25m**

Catalog Code: DEX-DET-B-025-025-25

Conditions:

Injection: Split 250°C, 0.5 μ L, 1:50 Split Ratio.

Sample dilution: 1mg/1mL each in n-Hexane.

Oven Temp.: 75°C, 1°C/min, 160°C.

Detector: FID 250°C.

Carrier Gas: Hydrogen, 75 kPa.

Peaks:

1. (+/-) Alpha-Pinene
2. (+/-) Linalool
3. (+/-) Menthol

