



Column: **MEGA-WAX - 0.32mm, 1.00 $\mu$ m, 60m**  
Catalog Code: S-WAX-032-100-60

**Conditions:**

Injection: Split 250°C, 0.5  $\mu$ L injection volume, 250 mL/min split flow.  
Detector: FID 250°C.  
Oven Program: 40°C, 8°C/min, 200°C.  
Carrier Gas: Helium, 150 kPa.

1. Hexane
2. Iso-Octane
3. Acetone
4. Methyl acetate
5. Ethyl acetate
6. 2-Butanone (MEK)
7. 2-Propanol (Isopropanol)
8. Ethanol
9. Trichloroethylene
10. Isobutyl acetate
11. Toluene
12. n-Butyl acetate
13. Isobutanol
14. 1-Methoxy-2-Propanol
15. Ethylbenzene
16. m-Xylene
17. p-Xylene
18. Isopropylbenzene
19. Styrene
20. 2-Ethoxyethyl acetate
21. Cyclohexanone
22. 2-Butoxyethanol
23. 2-Butoxyethyl acetate



since  
1980



contact us: [info@mega.mi.it](mailto:info@mega.mi.it)

 follow us [@MEGAColumns](https://twitter.com/MEGAColumns)

[www.mega.mi.it](http://www.mega.mi.it)